

Apolipoprotein A-IV

Catalog No: #AP80401



Package Size: #AP80401-1 50ug #AP80401-2 100ug #AP80401-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein A-IV
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Apo-AIV, Apolipoprotein A4
Accession No.	Uniprot ID:Q32PJ2
Uniprot	Q32PJ2
GeneID	537301;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APOA4

Note: This product is for in vitro research use only