

Aquaporin-1

Catalog No: #AP80402



Package Size: #AP80402-1 50ug #AP80402-2 100ug #AP80402-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aquaporin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AQP-1,Aquaporin-CHIP,Delayed early response protein 2,DER2,Water channel protein for red blood cells and kidney proximal tubule
Accession No.	Uniprot ID:Q02013
Uniprot	Q02013
GeneID	11826;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aqp1

Note: This product is for in vitro research use only