

## Alpha-hemoglobin-stabilizing protein

Catalog No: #AP80409



Package Size: #AP80409-1 50ug #AP80409-2 100ug #AP80409-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-hemoglobin-stabilizing protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Erythroid differentiation-related factor,Erythroid-associated factor,EDRF,ERAF
Accession No.	Uniprot ID:Q9NZD4
Uniprot	Q9NZD4
GeneID	51327;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AHSP

Note: This product is for in vitro research use only