

## A-kinase anchor protein 1, mitochondrial

Catalog No: #AP80410



Package Size: #AP80410-1 50ug #AP80410-2 100ug #AP80410-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A-kinase anchor protein 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	A-kinase anchor protein 149 kDa,AKAP 149,Dual specificity A-kinase-anchoring protein 1,D-AKAP-1,Protein kinase A-anchoring protein 1,PRKA1,Spermatid A-kinase anchor protein 84,S-AKAP84,AKAP149,PRKA1
Accession No.	Uniprot ID:Q92667
Uniprot	Q92667
GeneID	8165;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AKAP1

Note: This product is for in vitro research use only