

E3 ubiquitin-protein ligase AMFR

Catalog No: #AP80412



Package Size: #AP80412-1 50ug #AP80412-2 100ug #AP80412-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase AMFR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Autocrine motility factor receptor,AMF receptor,RING-type E3 ubiquitin transferase AMFR
Accession No.	Uniprot ID:Q9R049
Uniprot	Q9R049
GeneID	23802;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Amfr

Note: This product is for in vitro research use only