

Natriuretic peptides A

Catalog No: #AP80414



Package Size: #AP80414-1 50ug #AP80414-2 100ug #AP80414-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Cardiodilatin,CDD,Cardiodilatin-related peptide,CDP,Prepronatriodilatin
Accession No.	Uniprot ID:P24259
Uniprot	P24259
GeneID	397496;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NPPA

Note: This product is for in vitro research use only