

Apolipoprotein C-II

Catalog No: #AP80416



Package Size: #AP80416-1 50ug #AP80416-2 100ug #AP80416-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Apolipoprotein C-II |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Apo-CII, Apolipoprotein C2 |
| Accession No. | Uniprot ID:P19034 |
| Uniprot | P19034 |
| GeneID | 618039; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APOC2

Note: This product is for in vitro research use only