

Serum albumin

Catalog No: #AP80419



Package Size: #AP80419-1 50ug #AP80419-2 100ug #AP80419-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serum albumin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Alpha-livetin
Accession No.	Uniprot ID:P19121
Uniprot	P19121
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALB

Note: This product is for in vitro research use only