

Protein ALEX

Catalog No: #AP80423



Package Size: #AP80423-1 50ug #AP80423-2 100ug #AP80423-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein ALEX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alternative gene product encoded by XL-exon,Gnas1
Accession No.	Uniprot ID:Q9Z213
Uniprot	Q9Z213
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gnas

Note: This product is for in vitro research use only