

Ankyrin repeat domain-containing protein 30A

Catalog No: #AP80425

Package Size: #AP80425-1 50ug #AP80425-2 100ug #AP80425-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ankyrin repeat domain-containing protein 30A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Serologically defined breast cancer antigen NY-BR-1
Accession No.	Uniprot ID:Q9BXX3
Uniprot	Q9BXX3
GeneID	91074;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANKRD30A

Note: This product is for in vitro research use only