

# Annexin A6

Catalog No: #AP80426



Package Size: #AP80426-1 50ug #AP80426-2 100ug #AP80426-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Annexin A6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Annexin VI,Annexin-6,Calcium-binding protein 65/67,CBP 65/67,Anx6
Accession No.	Uniprot ID:P48037
Uniprot	P48037
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Anxa6

Note: This product is for in vitro research use only