

Natriuretic peptides B

Catalog No: #AP80427



Package Size: #AP80427-1 50ug #AP80427-2 100ug #AP80427-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Gamma-brain natriuretic peptide
Accession No.	Uniprot ID:P16859
Uniprot	P16859
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NPPB

Note: This product is for in vitro research use only