

# Annexin A1

Catalog No: #AP80428



Package Size: #AP80428-1 50ug #AP80428-2 100ug #AP80428-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Annexin A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Annexin I,Annexin-1,Calpactin II,Calpactin-2,Chromobindin-9,Lipocortin I,Phospholipase A2 inhibitory protein,p35,ANX1
Accession No.	Uniprot ID:Q92108
Uniprot	Q92108
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ANXA1

Note: This product is for in vitro research use only