

AKT-interacting protein

Catalog No: #AP80430



Package Size: #AP80430-1 50ug #AP80430-2 100ug #AP80430-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AKT-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Fused toes protein homolog,Fts
Accession No.	Uniprot ID:Q5FVH4
Uniprot	Q5FVH4
GeneID	291906;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aktip

Note: This product is for in vitro research use only