

Angiogenin

Catalog No: #AP80432



Package Size: #AP80432-1 50ug #AP80432-2 100ug #AP80432-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiogenin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiogenin-1,Ribonuclease 5,RNase 5,Ang1,Rnase5,Rnase5a
Accession No.	Uniprot ID:P21570
Uniprot	P21570
GeneID	11727;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ang

Note: This product is for in vitro research use only