

Adenylate cyclase type 10

Catalog No: #AP80438



Package Size: #AP80438-1 50ug #AP80438-2 100ug #AP80438-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AH-related protein,Adenylate cyclase homolog,Germ cell soluble adenylyl cyclase,hsAC,Testicular soluble adenylyl cyclase,SAC
Accession No.	Uniprot ID:Q96PN6
Uniprot	Q96PN6
GeneID	55811;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADCY10

Note: This product is for in vitro research use only