

# Beta-2 adrenergic receptor

Catalog No: #AP80439



Package Size: #AP80439-1 50ug #AP80439-2 100ug #AP80439-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-2 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-2 adrenoreceptor,Beta-2 adrenoceptor,Adrb2r
Accession No.	Uniprot ID:P10608
Uniprot	P10608
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Adrb2

Note: This product is for in vitro research use only