

α -tubulin (Acetyl Lys40) Monoclonal Antibody(4A8)(4A8)

Catalog No: #HW203

Package Size: #HW203-1 50ul #HW203-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	α -tubulin (Acetyl Lys40) Monoclonal Antibody(4A8)(4A8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4A8
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	α -tubulin (Acetyl Lys40) Mouse Monoclonal Antibody(4A8) detects endogenous levels of α -tubulin (Acetyl Lys40)
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	α -tubulin
Accession No.	Swiss-Prot:P68363Gene ID:10376
Uniprot	P68363
GeneID	10376;
Calculated MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:50-1:100


Images



Western blot analysis of extracts from HeLa cells, untreated (-) or treated with TSA (1 μ M, 18 hr; +), using Acetyl- α -tubulin(Lys40) Mouse mAb 1:2000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using α -tubulin(Acetyl Lys40) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using α -tubulin(Acetyl Lys40) Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only