

## Aldehyde dehydrogenase, cytosolic 1

Catalog No: #AP80440



Package Size: #AP80440-1 50ug #AP80440-2 100ug #AP80440-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aldehyde dehydrogenase, cytosolic 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ALDH class 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A7,Aldehyde dehydrogenase phenobarbital-inducible,Aldh-pb
Accession No.	Uniprot ID:O35945
Uniprot	O35945
GeneID	26358;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Aldh1a7

Note: This product is for in vitro research use only