

Natriuretic peptides B

Catalog No: #AP80441



Package Size: #AP80441-1 50ug #AP80441-2 100ug #AP80441-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Gamma-brain natriuretic peptide
Accession No.	Uniprot ID:P13204
Uniprot	P13204
GeneID	508734;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NPPB

Note: This product is for in vitro research use only