

Androgen receptor

Catalog No: #AP80446



Package Size: #AP80446-1 50ug #AP80446-2 100ug #AP80446-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Androgen receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dihydrotestosterone receptor,Nuclear receptor subfamily 3 group C member 4,Nr3c4
Accession No.	Uniprot ID:P19091
Uniprot	P19091
GeneID	11835;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ar

Note: This product is for in vitro research use only