

AMPK α 1 Monoclonal Antibody(5G11)

Catalog No: #44102

Package Size: #44102-1 50ul #44102-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	AMPK α 1 Monoclonal Antibody(5G11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5G11
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	AMPK α 1 Mouse Monoclonal Antibody(5G11) detects endogenous levels of AMPK α 1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	AMPK α 1
Accession No.	Swiss-Prot:Q13131Gene ID:5562
Uniprot	Q13131
GeneID	5562;
Calculated MW	63kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:50-1:100

Images

Western blot analysis of 1)Hela, 2) 293T, 3)3T3, 4) PC12 with AMPK α 1 Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse ColonTissue using AMPK a1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only