

## Membrane primary amine oxidase

Catalog No: #AP80451



Package Size: #AP80451-1 50ug #AP80451-2 100ug #AP80451-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Membrane primary amine oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Copper amine oxidase, Semicarbazide-sensitive amine oxidase, SSAO, Vascular adhesion protein 1, VAP-1
Accession No.	Uniprot ID:Q9TTK6
Uniprot	Q9TTK6
GeneID	281002;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: AOC3

Note: This product is for in vitro research use only