

AP-1 complex subunit mu-1

Catalog No: #AP80454

Package Size: #AP80454-1 50ug #AP80454-2 100ug #AP80454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AP-1 complex subunit mu-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AP-mu chain family member mu1A,Adaptor protein complex AP-1 subunit mu-1,Adaptor-related protein complex 1 subunit mu-1,Clathrin assembly protein complex 1 mu-1 medium chain 1,Golgi adaptor HA1/AP1 adaptin mu-1 subunit,Mu-adaptin 1,Mu1A-adaptin
Accession No.	Uniprot ID:Q2KJ81
Uniprot	Q2KJ81
GeneID	504310;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AP1M1

Note: This product is for in vitro research use only