

# Jouberin

Catalog No: #AP80456



Package Size: #AP80456-1 50ug #AP80456-2 100ug #AP80456-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Jouberin  |
| Brief Description  | Recombinant Protein                                     |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                            |
| Species Reactivity | Rt  |
| Other Names        | Abelson helper integration site 1 protein homolog,AHI-1 |
| Accession No.      | Uniprot ID:Q6DTM3                                       |
| Uniprot            | Q6DTM3  |
| GeneID             | 308923;   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ahi1

Note: This product is for in vitro research use only