

## Single-stranded DNA cytosine deaminase

Catalog No: #AP80457



Package Size: #AP80457-1 50ug #AP80457-2 100ug #AP80457-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Single-stranded DNA cytosine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase
Accession No.	Uniprot ID:Q2PT36
Uniprot	Q2PT36
GeneID	539888;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AICDA

Note: This product is for in vitro research use only