

Interferon-inducible protein AIM2

Catalog No: #AP80459



Package Size: #AP80459-1 50ug #AP80459-2 100ug #AP80459-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Interferon-inducible protein AIM2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Absent in melanoma 2 |
| Accession No. | Uniprot ID:O14862 |
| Uniprot | O14862 |
| GeneID | 9447; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AIM2

Note: This product is for in vitro research use only