

## AKT-interacting protein

Catalog No: #AP80463



Package Size: #AP80463-1 50ug #AP80463-2 100ug #AP80463-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AKT-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FT1,Fused toes protein,Fif,Ft1,Fts
Accession No.	Uniprot ID:Q64362
Uniprot	Q64362
GeneID	14339;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Aktip

Note: This product is for in vitro research use only