

Aminopeptidase O

Catalog No: #AP80464



Package Size: #AP80464-1 50ug #AP80464-2 100ug #AP80464-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aminopeptidase O
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AP-O, Onpep
Accession No.	Uniprot ID:Q8BXQ6
Uniprot	Q8BXQ6
GeneID	72061;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Aoep

Note: This product is for in vitro research use only