

Type-2 angiotensin II receptor

Catalog No: #AP80468



Package Size: #AP80468-1 50ug #AP80468-2 100ug #AP80468-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Type-2 angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin II type-2 receptor,AT2
Accession No.	Uniprot ID:P35374
Uniprot	P35374
GeneID	11609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Agtr2

Note: This product is for in vitro research use only