

Allograft inflammatory factor 1

Catalog No: #AP80469



Package Size: #AP80469-1 50ug #AP80469-2 100ug #AP80469-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Allograft inflammatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AIF-1,Ionized calcium-binding adapter molecule 1,MRF-1,Microglia response factor,Iba1,Mrf1
Accession No.	Uniprot ID:P55009
Uniprot	P55009
GeneID	29427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aif1

Note: This product is for in vitro research use only