

I κ B β Monoclonal Antibody(8D11)

Catalog No: #44104

Package Size: #44104-1 50ul #44104-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	I κ B β Monoclonal Antibody(8D11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8D11
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	I κ B β Mouse Monoclonal Antibody(8D11) detects endogenous levels of I κ B β
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	I κ B β
Accession No.	Swiss-Prot:Q15653Gene ID:4793
Uniprot	Q15653
GeneID	4793;
Calculated MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using I κ B β (Mouse mAb diluted at 1:200).



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using $\text{I}\kappa\text{B}$ β Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only