

Aldose reductase

Catalog No: #AP80470



Package Size: #AP80470-1 50ug #AP80470-2 100ug #AP80470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldose reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AR,Aldehyde reductase,Akr1b4,Aldr1
Accession No.	Uniprot ID:P07943
Uniprot	P07943
GeneID	24192;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akr1b1

Note: This product is for in vitro research use only