

Proline-rich AKT1 substrate 1

Catalog No: #AP80471



Package Size: #AP80471-1 50ug #AP80471-2 100ug #AP80471-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proline-rich AKT1 substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proline-rich AKT substrate,Pras
Accession No.	Uniprot ID:Q9D1F4
Uniprot	Q9D1F4
GeneID	67605;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Akt1s1

Note: This product is for in vitro research use only