

Aldehyde dehydrogenase, cytosolic 1

Catalog No: #AP80472



Package Size: #AP80472-1 50ug #AP80472-2 100ug #AP80472-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase, cytosolic 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ALDH class 1,ALDH-E1,ALHDII,Aldehyde dehydrogenase family 1 member A7,Aldehyde dehydrogenase phenobarbital-inducible,Aldh-pb,Aldh1,Aldh1a4
Accession No.	Uniprot ID:P13601
Uniprot	P13601
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldh1a7

Note: This product is for in vitro research use only