

Alpha-actinin-1

Catalog No: #AP80474



Package Size: #AP80474-1 50ug #AP80474-2 100ug #AP80474-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-actinin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-actinin cytoskeletal isoform,F-actin cross-linking protein,Non-muscle alpha-actinin-1
Accession No.	Uniprot ID:Q9Z1P2
Uniprot	Q9Z1P2
GeneID	81634;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Actn1

Note: This product is for in vitro research use only