

Anti-Muellerian hormone type-2 receptor

Catalog No: #AP80476



Package Size: #AP80476-1 50ug #AP80476-2 100ug #AP80476-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Muellerian hormone type-2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Anti-Muellerian hormone type II receptor,AMH type II receptor,MIS type II receptor,MISRII,AMHR,MISR2
Accession No.	Uniprot ID:Q16671
Uniprot	Q16671
GeneID	269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AMHR2

Note: This product is for in vitro research use only