

Actin, aortic smooth muscle

Catalog No: #AP80479



Package Size: #AP80479-1 50ug #AP80479-2 100ug #AP80479-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Actin, aortic smooth muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-actin-2,Actsa,Actvs
Accession No.	Uniprot ID:P62737
Uniprot	P62737
GeneID	11475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acta2

Note: This product is for in vitro research use only