

## Aldehyde dehydrogenase, dimeric NADP-preferring

Catalog No: #AP80484



Package Size: #AP80484-1 50ug #AP80484-2 100ug #AP80484-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aldehyde dehydrogenase, dimeric NADP-preferring
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Aldehyde dehydrogenase family 3 member A1,HTC-ALDH,Tumor-associated aldehyde dehydrogenase,Aldd,Aldh3
Accession No.	Uniprot ID:P11883
Uniprot	P11883
GeneID	25375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Aldh3a1

Note: This product is for in vitro research use only