

Amelogenin, X isoform

Catalog No: #AP80486



Package Size: #AP80486-1 50ug #AP80486-2 100ug #AP80486-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amelogenin, X isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Leucine-rich amelogenin peptide,LRAP,Amel,Amgx
Accession No.	Uniprot ID:P63278
Uniprot	P63278
GeneID	29160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Amelx

Note: This product is for in vitro research use only