

# Methionine aminopeptidase 2

Catalog No: #AP80488



Package Size: #AP80488-1 50ug #AP80488-2 100ug #AP80488-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Methionine aminopeptidase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP 2, MetAP 2, Initiation factor 2-associated 67 kDa glycoprotein, Peptidase M 2, p67eIF2, p67, MNPEP, P67EIF2
Accession No.	Uniprot ID: P50579
Uniprot	P50579
GeneID	10988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: METAP2

Note: This product is for in vitro research use only