

Angiopoietin-1

Catalog No: #AP80489



Package Size: #AP80489-1 50ug #AP80489-2 100ug #AP80489-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ANG-1,Agpt
Accession No.	Uniprot ID:O08538
Uniprot	O08538
GeneID	11600;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Angpt1

Note: This product is for in vitro research use only