

Angiopoietin-2

Catalog No: #AP80490



Package Size: #AP80490-1 50ug #AP80490-2 100ug #AP80490-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ANG-2
Accession No.	Uniprot ID:A0A8J8
Uniprot	A0A8J8
GeneID	607616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANGPT2

Note: This product is for in vitro research use only