

Angiotensin-like protein 1

Catalog No: #AP80496



Package Size: #AP80496-1 50ug #AP80496-2 100ug #AP80496-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiotensin-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q8IY63
Uniprot	Q8IY63
GeneID	154810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AMOTL1

Note: This product is for in vitro research use only