

Aldo-keto reductase family 1 member C2

Catalog No: #AP80498

Package Size: #AP80498-1 50ug #AP80498-2 100ug #AP80498-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Aldo-keto reductase family 1 member C2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | 3-alpha-HSD3, Chlordecone reductase homolog HAKRD, Dihydrodiol dehydrogenase 2, DD-2, Dihydrodiol dehydrogenase/bile acid-binding protein, DD/BABP, Trans-1,2-dihydrobenzene-1,2-diol dehydrogenase, Type III 3-alpha-hydroxysteroid dehydrogenase, DDH2 |
| Accession No. | Uniprot ID:P52895 |
| Uniprot | P52895 |
| GeneID | 1646; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AKR1C2

Note: This product is for in vitro research use only