

Alpha-1A adrenergic receptor

Catalog No: #AP80500



Package Size: #AP80500-1 50ug #AP80500-2 100ug #AP80500-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-1A adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-1A adrenoreceptor,Alpha-1A adrenoceptor,Alpha-1C adrenergic receptor,Adra1c
Accession No.	Uniprot ID:P97718
Uniprot	P97718
GeneID	11549;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adra1a

Note: This product is for in vitro research use only