

## Adiponectin

Catalog No: #AP80501



Package Size: #AP80501-1 50ug #AP80501-2 100ug #AP80501-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Adiponectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	30 kDa adipocyte complement-related protein,Adipocyte complement-related 30 kDa protein,ACRP30,Adipocyte, C1q and collagen domain-containing protein,Adipose most abundant gene transcript 1 protein,apM-1,ACRP30,APM1
Accession No.	Uniprot ID:Q3Y5Z3
Uniprot	Q3Y5Z3
GeneID	282865;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADIPOQ

Note: This product is for in vitro research use only