

Protein argonaute-2

Catalog No: #AP80502



Package Size: #AP80502-1 50ug #AP80502-2 100ug #AP80502-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Protein argonaute-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Argonaute2, Argonaute RISC catalytic component 2, Eukaryotic translation initiation factor 2C 2, eIF-2C 2, Golgi ER protein 95 kDa, GERp95, Protein slicer, Eif2c2 |
| Accession No. | Uniprot ID:Q9QZ81 |
| Uniprot | Q9QZ81 |
| GeneID | 59117; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ago2

Note: This product is for in vitro research use only