

Atrial natriuretic peptide receptor 1

Catalog No: #AP80503



Package Size: #AP80503-1 50ug #AP80503-2 100ug #AP80503-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Atrial natriuretic peptide receptor type A, ANP-A, Guanylate cyclase A, GC-A
Accession No.	Uniprot ID:P18910
Uniprot	P18910
GeneID	24603;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Npr1

Note: This product is for in vitro research use only