

## DNA-(apurinic or apyrimidinic site) lyase

Catalog No: #AP80505



Package Size: #AP80505-1 50ug #AP80505-2 100ug #AP80505-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA-(apurinic or apyrimidinic site) lyase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APEX nuclease, APEN, Apurinic-apyrimidinic endonuclease 1, AP endonuclease 1, REF-1, Redox factor-1, APE, APE1, APEX, APX, HAP1, REF1
Accession No.	Uniprot ID:P27695
Uniprot	P27695
GeneID	328;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: APEX1

Note: This product is for in vitro research use only